



North-Eastern Hill University
Department of Chemistry
NEHU Campus, Shillong, 793022 (Meghalaya)

Rishanlang Nongkhlaw, Ph.D, Assistant Professor

Tel: 0364-272-2628 (O)/2631(Lab); 9436101349(M)

E-mail: rlnongkhlaw@nehu.ac.in/rlnongkhlaw@gmail.com

Ref No: CHEM/RLN/SERB/2017/002

Date: 24.03.2017

NOTICE INVITING TENDER

Sealed quotations are invited from reputed manufacturers/ authorized dealers for supply of instruments for the DST-SERB funded project in the Department of Chemistry, North-Eastern Hill University, Shillong. The instruments along with specifications are listed in the **Annexure-I**. The quotation should reach the undersigned by post on or before the **15th April, 2017 (up to 5:00 pm)**. Quotation received after the due date and time will not be entertained.

The undersigned reserves the right to accept or reject the tender without stating any reason whatsoever.

Terms and Conditions:

1. The supply should be completed within 30 days from the date of issue of supply.
2. The payment will be completed only after the installation of the instruments.
3. The instruments should be delivered and installed without any additional cost.
4. Users list should be attached along with the quotation.
5. The company should provide factory trained engineers for maintenance.

Dr. Rishanlang Nongkhlaw,
DST-SERB funded project, Principal
Investigator, Department of Chemistry,
North-Eastern Hill University
Mawkynroh, Umshing, Shillong - 793022
Email: rlnongkhlaw@gmail.com;
rlnongkhlaw@nehu.ac.in

ANNEXURE-1

Sl. No.	Instruments	Specifications	Quantity
1	Electronic Balance	Sartorius Electronic Balance: <ul style="list-style-type: none">• Model Quintix224-10IN.• Capacity of 220 g and readability of 0.1 mg.• Pan size 90 mm diameter.• Comes with touch screen and complete instruction manual.	1
2	Centrifuge	Centrifuge Remi Brushless Motor R8M BL 16 tubes: <ul style="list-style-type: none">• Capacity of 16 tubes.• Stable speed output even under unstable voltage conditions.• Brushless induction motor with variable frequency drive and microprocessor controller and digital display.• Inverter fault detection system with auto shutdown.• Safety lid interlock to prevent lid opening during centrifugation along with imbalance detection & centrifugation stop with display of error.	1
3	Ultrasound Instrument	Sonics VC 505: <ul style="list-style-type: none">• 500 Watt Ultrasonic Processor.• 250 microliters to 1 liter.	1
4	Chiller and Easy Vac vacuum control	Chiller and Easy Vac vacuum control: <ul style="list-style-type: none">• Chiller for Rotavapor (Temperature range -25°C to 80°C)• Easy-Vac for Buchi vacuum pump.	1

Quotation should be posted to:

Dr. Rishanlang Nongkhaw,
DST-SERB funded project, Principal Investigator,
Department of Chemistry,
North-Eastern Hill University Mawkynroh,
Umshing, Shillong - 793022
Email: rlnongkhaw@gmail.com;
rlnongkhaw@nehu.ac.in